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GB How to clean your teeth properly

Join the first animated, practical instruction according to the "New Swiss School of Oral hygiene" for a correct use of toothbrush, floss and interproximal brushes. Visit the www.curaden.ch homepage to learn how to maintain fresh breath, healthy gums and teeth without plaque and stain.

1. Hold your toothbrush clearly angulated against the gums. One half of the bristles covers the gums, the other half covers the surface of the teeth.
2. The brush is slowly guided along the tooth-row from the front to the back by applying little pressure and fine circular movements.
3. After the careful cleaning of the outside of the first half of the upper jaw you proceed to the other half.
4. Proceed to cleaning the lingual surface of the upper teeth as for the buccal surfaces of the lower teeth.
5. After the buccal cleansing, lingual cleansing is on.
6. The compact head of the CS 5460 ideally suits the vertical access to difficult to reach lingual anterior areas.
7. While brushing the rear surfaces of the last molars the advantage of the compact head becomes particularly evident.
8. The compact head allows perfect lingual cleansing even in case of narrow jaw anatomy.
9. Last but not least. Brush all occlusal surfaces, particularly important for children.

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